



# MERCER I

## 60" x 32" Drop-In Acrylic Whirlpool



Model Number  
**3260MW064**  
Drop-In Whirlpool in White

### Features



- Color-matched jet trim
- 10 adjustable hydrotherapy jets
- Pillow included
- Constructed of high-gloss durable acrylic
- 3-Speed control for the 9.0 Amp pump
- Inline heater available – LM300
- Mood light available – LM522
- Pre-leveled base for easy installation
- Made in the USA



**\*All dimensions are +/- 1/2" and subject to change without notice\***  
**Confirm specifications before installation**



Standard (S) Pump Location and  
Alternate (A-1 & A-2) Pump locations.



### Technical Information

- Cut-Out Dimensions: 57 3/4" x 29 7/8"
- Weight: 120 lbs.
- Water Depth to Overflow: 15"
- Water Operating Level: 44 gal.
- Water Capacity: 54 gal.
- Sump Dimension: 37" x 16"

### Electrical Requirements.

- Standard Pump: 9.0 Amps
- (Optional) LM501 Pump Upgrade: 13.5 Amps
- (Optional) LM300 Heater: 12.0 Amps
- (Optional) LM509 Back Massage: 7.5 Amps
- Dedicated GFCI's required for each option

# MERCER I

**3260MW064**  
Drop-In Whirlpool in White

## Whirlpool

Experience a powerful, deep tissue therapeutic massage, which will soothe muscles & joints while releasing tension. Our eight standard, high performance jets will promote increased blood flow, improved flexibility and help reduce accumulated toxins.

Laurel Mountain enables you to enhance your bathing experience by offering a broad range of options to customize your bath. The options available on this bath are shown below.

## AVAILABLE OPTIONS & UPGRADES

For Pricing, visit the “Help” page at [www.laurelmountainwhirlpools.com](http://www.laurelmountainwhirlpools.com). Download the Price book. For pictures and descriptions of options, go to [www.laurelmountainwhirlpools.com](http://www.laurelmountainwhirlpools.com) and click “options” tab found at the bottom of the website.

### Therapeutic

<b>LM300 Inline Heater</b>	<b>BACK MASSAGE SYSTEM</b>	LM509 White
<b>LM522 Chromatherapy</b>	LM505 Chrome	LM510 Biscuit
<b>LM500 Extra Jets (x2)</b>	LM506 Brushed Nickel	LM511 Almond
<b>LM515 Neck Massage Pillow</b>	LM507 Bronze	LM512 Linen
	LM508 Polished Brass	LM513 Bone

### Accessories

<b>TRIM KITS</b>	<b>DEEP SOAK DRAINS</b>	<b>8” GRAB BARS</b>
LM101 Chrome	LM200 Chrome	LM230 White
LM102 Brushed Nickel	LM201 Brushed Nickel	LM231 Biscuit
LM103 Bronze	LM202 Bronze	LM232 Almond
LM104 Polished Brass	LM203 Polished Brass	LM233 Linen
	LM210 White	LM234 Bone
<b>LM321 Electronic Control, Timer, &amp; Variable Speed Pump</b>	LM211 Biscuit	<b>LM335 Ozonator (Sanitizing System)</b>
	LM212 Almond	
	LM213 Linen	
	LM214 Bone	

### Design / Construction

<b>LM537 Bariatric/Heavy Duty</b>	<b>LM542 Plastic Flange</b>	<b>LM555 Skirt (RH drain)</b>
<b>LM540 Tile Flange (RH)</b>	<b>LM545 A-1 Pump Loc</b>	<b>LM556 Skirt (LH drain)</b>
<b>LM541 Tile Flange (LH)</b>	<b>LM546 A-2 Pump Loc</b>	

## INSTALLATION

Install this unit according to the Owner’s Manual. The Owner’s manual can be found by visiting the “Help” page at [www.laurelmountainwhirlpools.com](http://www.laurelmountainwhirlpools.com). Under “Relevant Links” on the left side of the page, click “Owner’s Manual”. Then click “Colony Owner’s Manual”.